

HOME BEAUTY INSTRUMENT EXPERT

PRODUCT MANUAL

HOME BEAUTY INSTRUMENT EXPERT

Shenzhen Kizoku Life Technology Co., Ltd.

Add : 3F South, Building B, No.71-2 Xintian Road, Fuyong Street,
Bao'an District, Shenzhen, 518103 China.

Tel : +86-755-2373-2451

Fax : +86-755-2333-8247

Web : www.nt-beauty.com

APP QR



Web QR



Wechat QR





We can't change the time
but we try to mitigate the damage to your beauty with
future technology and products.



CONTENTS



01

Beauty and
Skin Care
01–16



02

Skin Cleansing
Care
17–22



03

Hydrate and
Activate Skin
23–30




04

Body Shaping
31–34



05

Hair care &
styling
35–36



06

Healthy Life
37–44



07


New Arrivals
45–58



Smart Photon RF
Tender Skin
Appratus
·F81E



RF Rejuvenation, Photon
Beauty New Experience
of Intelligent Beauty Care



02. Multi-function synchronization effect of patent technology

- **RF+Red-ray**
Double skin renewal energy, enliven collagen vitality
- **RF+Sonic+ION exporting**
RF and Sonic energy open pores, ION adsorption realizes deep cleaning
- **RF+ENI+Red-ray+Sonic importing**
RF and Sonic energy open pores, pulse waveguide input assistance, photon skin renewal promotes absorption
- **RF+EMS+Red-ray+Blue-ray+Sonic lifting**
RF and Photon energy rejuvenate skin, sonic massage plays a role of physical beating, micro-current lifting and shaping
- **COOL+Blue-ray**
Ice wave shrinks pore, locks in moisture, blue ray reduces inflammation and calms skin

03. APP connection, new experience of family beauty care

NanoSkin has built-in Bluetooth, data storage, intelligent connection "Beauty Time" app, monitoring your beauty status at any time.



04. Multiple intelligent modes, easier to use

This product is specially equipped with four intelligent modes: photon beauty, photon nutrition, photon V-face and photon importing, so that users can operate more simply and intuitively.

Step	Function	Photon-Beauty	Photon-Nutrition	Photon-V-face	Photon-Import
Deep Cleansing	RF+ION Cotton pad +Toner	✓	✓	✓	✓
Import Assistant	RF+ENI RED-RAY Essence/ mask	✓	✓	✓	✓
V-face	RF+EMS PHOTON Gel/moisturizing cream	✓	✓	✓	✓
Cold Compress	COOL BLUE-RAY Direct use	✓	✓	✓	✓

✦ SPECIFICATIONS

Net weight:	About 300g	Battery Specification:	7.4V, 650mAh	Rated Power:	About 10W
L C D:	1.5 inches, 240*240	Material:	ABS, Aluminum Alloy, Stainless Steel	Full-charge Time:	About 2 hours
Dimensions:	H190×W55×D64.5mm	Adapter:	Input 100-240V, 0.35A Output 9V=1.5A		

Photon RF Tender
Skin Apparatus

•F82E



Radio Frequency
Rejuvenation
Photon Beauty



01.

Integrating various cosmetic techniques originating from professional fields

RF technology from the professional field with ION export, ENI import, EMS facial lifting, photon rejuvenation, sonic massage, COOL compress, so that you can enjoy the new experience of professional beauty care at home.

02. Multi-function synchronization effect of patent technology

- **RF+Red-ray**
Double skin renewal energy, enliven collagen vitality
- **RF+Sonic+ION exporting**
RF and Sonic energy open pores, ION adsorption realizes deep cleaning
- **RF+ENI+Red-ray+Sonic importing**
RF and Sonic energy open pores, pulse waveguide input assistance, photon skin renewal promotes absorption
- **RF+EMS+Red-ray+Blue-ray+Sonic lifting**
RF and Photon energy rejuvenate skin, sonic massage plays a role of physical beating, micro-current lifting and shaping
- **COOL+Blue-ray**
Ice wave shrinks pore, locks in moisture, blue-ray reduces inflammation and calms skin



✦ SPECIFICATIONS

Net weight: About 300g	Battery Specification: 7.4V, 650mAh	Rated power: About 10W
L C D: 1.5 inches	Material: ABS, Aluminum Alloy, Stainless Steel	Full-charge time: About 2 hours
Dimensions: H190×W55×L64.5mm	Adapter: Input 100-240V, 0.35A	Output 9V=1.5A

nanoSkin-L nanoSkin-D
RF beauty technique
Ultimate cool experience



Luxury RF tender
skin apparatus.

•F82



luxury RF tender
skin apparatus
(cool version)

•F82C

01. Innovative integration, multiple functions for one device



RF-ION
cleansing and
detoxifying



RF-ENI
Deep
absorption




RF-EMS
Face-lifting
and shaping




RF
Skin
rejuvenation



02. Professional skin care of beauty salon level, multiple effects for your beauty care



Lack of
collagen
and elastin



Refresh
collagen
and elastin



03. Advanced double beauty functions working at a time, better effects for your beauty care

The dermis would be warmed up by the benefit of RF function, which can improve the activity and metabolism of collagen cells and open the pores. Meanwhile, the functions such as ION deep cleaning, ENI electroporation importing, EMS electrical muscle stimulation would be applied together as three smart combined mode for Double Beauty function to intensify the effects of beauty care.



- **RF:** RF beauty function can activate collagen, makes your skin smooth and tender, also improves dullness and helps you anti-aging.
- **EMS:** The micro-current muscle stimulation technology can help in exercising deep facial muscles, lifting and tightening your face.
- **ION:** ION exporting, thoroughly clean away the dirt and make-up residue from inside to outside.
- **ENI:** Electroporation Importing, penetrate nutrients straight into the cells for deep and highly-efficient absorption.
- **COOL:** Tightening the pores, tighten your skin by cold compress, also help retain moisture and nutrients.

Note: F82C has this functions, F82 is not.

⚡ SPECIFICATIONS

Battery voltage: 7.4V =	Battery capacity: 650mAh	Screen size: 1.5 inches
Input: 9V=1.5A	Net weight: About 280g	Gear: 1-5 adjustable
Rated power: About 10W	Dimensions: H190×W55×D64.5mm	Full-charge time: About 2 hours
Material: ABS, Aluminum Alloy, Stainless Steel		

nanoskin
Smart RF tender
skin apparatus



nanoskin
•F81



01. Four intelligent modes, double beauty function synchronous function

Innovative RF+ION, RF+ENI, RF+EMS, dual function synchronous mode, mobile phone app connection, synchronous management, intelligent design, designed for simple operation of home beauty care.



Smart beauty mode
RF+ION
RF+ENI
RF+EMS
Cool



Nutritious mode
RF+ION
RF+ENI
Cool



Smart V face-lifting mode
RF+ EMS
Cool



Nutrition lock
RF+ ENI
Cool



02. Innovation integration One machine with 5 purposes



RF
Skin
rejuvenation



EMS
RF+EMS
Face-lifting and
shaping



ION
RF+ION
cleansing and
detoxifying



COOL
COOL
Ice compress
and tightening



ENI
RF+ENI
Deep
absorption



SPECIFICATIONS

Battery voltage: 7.4V =	Battery capacity: 650mAh	Screen size: 1.5 inches
Input: 9V=1.5A	Net weight: About 280g	Gear: 1-5 adjustable
Bluetooth: Support mobile app connection	Screen pixel: 240*240	Full-charge time: About 2 hours
Material: ABS, Aluminum Alloy, Stainless Steel	Dimensions: H190×W55×D64.5mm	Rated power: About 10W

**01. Innovative beauty function synchronization,
new experience and more functions eyes**sonic
massagephoton-
assistedhot compress
relaxation

ion importing


EMS
micro-current

ion exporting

trophoSkin-i innovatively integrates sonic massage, photon-assisted, hot compress relaxation, ion importing and exporting, EMS micro-current muscle exercise and other physical beauty care technology, according to the proposal of beauty experts for a scientific model combination, this brings a new beauty care experience.

**02. Multifunctional synchronous photon
importing mode to assist skin care products
absorption**

Ion importing function of trophoSkin-i is synchronous introduction with sonic massage, hot compress relaxes and opens pores to make skin care molecules easier to absorb, photon-assisted is synchronization, this brings a new beauty skin care experience.

**03. Multifunctional synchronous ion exporting
mode to clean pores and permeate skin**

TrophoSkin-i is unique with multidimensional exporting mode. In this mode, sonic massage and hot compress will open pores, ion function will absorb pore dirt to achieve deep cleaning, so that the skin is more transparent, alleviating acne problem, pore size and other issues.

**04. Micro-current and photon pulling mode
with multi function synchronous, shaping V
face.**

trophoSkin-i is built-in micro-current photon lifting mode. In this mode, EMS micro-current muscle exercise will lift expression muscle, hot compress function will activate the vitality of skin cells, photon relaxation function will work at the same time, four functions of synchronization will assist facial shaping, compact V face.

**✦ SPECIFICATIONS**

Battery Specification: 3.7V, 800mAh	Charger requirement: Output 5V=1A	Net weight: About 125g
Charging time: About 2-3 hours	Standby current: Less than 20 uA	Dimensions: H168 X W50 X L50mm
Consumed power: Maximally about 3 W	Material: ABS, PP, Stainless Steel	

Photon importing
instrument

-F65I



trophoSkin-i

Photon importing sonic
massage assisted pulling up
and making your skin comeliness

Ultrasonic importing
instrument

·F65U



Ultrasonic importing auxiliary absorption
one machine double effective

01. Innovative beauty synchronization function, leading market

Ultrasonic importing of trophoSkin-U will work with ion function synchronous, sonic massage, hot compress soothing and photon-assisted will bring a new beauty skin care experience, skin care nutrition are better absorption, but also to avoid the surface residual breeding bacteria bring other skin problems.



02. Ultrasonic pulling mode, shaping V face

trophoSkin-U is built-in ultrasonic lifting mode. In this mode, hot compress activates skin cells, ultrasonic and sonic massage physical lifting is synchronous, EMS micro-current muscle training is assistance, four synchronous functions will assist facial shaping and create delicate V-face.



03. Ion exporting , cleansing pores and permeating skin

trophoSkin-U is special of multi-dimensional ion exporting mode. In this mode, sonic massage and hot compress will open pores, ion exporting will absorb the dirt of pore to achieve deep cleaning, that makes the skin more permeable, so that follow-up care products is more easily to absorbed.



04. Photorejuvenation, a variety of beauty experience

- In **photon rejuvenation mode**, the red and blue light will alternate effect with sonic massage that will soothe the skin, and make it tender and light.
- In **photon massage mode**, hot compress and red light is synchronous effect, and the sonic massage will relax the skin
- In **photon relaxation mode**, hot compress and blue light is synchronize effect, and sonic massage will relieve the skin.



✦ SPECIFICATIONS

Battery Specification: 3.7V, 800mAh	Charger requirement: Output 5V=1A	Net weight: About 140g
Charging time: About 2-3 hours	Standby current: Less than 20 uA	Dimensions: H168×W50×L50mm
Consumed power: Maximally about 3 W	Material: ABS, PP, Stainless Steel	

Eyemagic

Magic Technology
focusing on
beautifying eyes

Eyemagic-E (lithium battery)


• F52E

Eyemagic-X (AA Battery)

• F52X

**01. Various cosmetic functions are affecting together ,that Eyemagic can be more effective**

Eyemagic uses the method innovatively that various cosmetic functions are affecting together, with scientific and technological magic to carry out eye care. It makes customers have better sense of body and better effects.

**02. Four nursing functions synchronous nutrition introduction Simultaneous massage and relaxation with three soothing functions**

40°C hot compress, sonic massage, blue light soothing, three technologies work together to massage your eyes and relax your skin around eyes.



40°C Constant Temperature



Sonic massage




Red light make skin tender




Ion exporting

Note: F52x has no 40°C hot compress

03. Nourishing and import mode**Ion-Import**


Ion-Import with constant temperature, and red light make skin tender. Sonic Massage function synchronization, greatly enhance the absorption of eye cream, improves the fat particle problem.

**Sonic massage**

Sonic massage relieves eye fatigue, promotes absorption of eye cream, and improves eye lines and dark under-eye circles.

**40°C Constant Temperature**

The 40°C constant temperature hot compress is specially designed for the skin of the eye. It is gentle and comfortable. It helps open the pores and improves the vitality and absorption of the cells in the eye skin.

**Red light make skin tender**

LED red light rejuvenates skin cells of the eye, and light energy supplements heat energy to further enhance the effect of the care.

**04. Smart design, easy to use**

- Smart application switch, starting when touching the skin, more energy saving;
- Smart switching mode, magic wand of beauty eye, easy to use.

Note: F52X has no such functions.

SPECIFICATIONS

Product Model: F52E	Hot temperature: about 40°C	Product Model: F52X
Net weight: about 43g	Charging time: about 1.5 hours	Net weight: about 32g
Material: ABS, PC	Built-in battery: 320mAh lithium battery	Dimensions: 144×23×21mm
Dimensions: 144×23×21mm	Charging method: 5V /1A or USB charging	Material: ABS, PC



V face
•F83

V face -L
•F83-L




Vface Vface-L

Roll the magic wand in your hand to build delicate V-shape face

01. Solar microcurrent, it can be used in the light without charging.

Solar panel on the handle can absorb microcurrent generated by the light ray and the microcurrent can be transmitted into skin by roller, promoting facial circulation to replenish human bioelectricity and restore skin elasticity.



02. 115°V-shape Gold Angle Massage Roller

Global surgery beauty experts found the beautiful mandibular angle is 115° after analyzing tens of thousands of female faces. Therefore, 115° is the prime angle for the most beautiful face in the surgery of mandibular.



03. Innovate four drainage massage roller that can fit each facial part from any angles.

It is equipped with exquisite balance interval and can generate four pressure points in sync by single rolling, achieving deeply seamless lifting for facial complex curve.

04. Be more beautiful with multiple purposes

Portable with small size, massage absorbing effect will be greater matching daily skin care product. Use together with mask, essence and cream. Enhancing the absorption effect for several times.

05. Deep lifting massage of SPA, reach women's natural beauty



Use for eye bag roll to temple from the skin below the eye in semi-circle way.



Use for crows feet roll in the small range from canthus to temples.



Vshape lifting scroll up repeatedly from mouth to ear to firm your face.



Vface lift Roll up repeatedly from chin near to ear.



Use for forehead gently roll from eyebrows to hairline.



Use for cervical Scroll up and down on both sides of neck and pay attention to avoid Adam's apple.

✦ SPECIFICATIONS

Product Model: F83L	Roller: 2 Rollers	Product Model: F83	Roller: 4 Rollers
Current: 15µA	Material: Aluminum Alloy, Refined Copper	Current: 15µA	Material: Aluminum Alloy, Refined Copper
Net weight: 118g	Dimensions: H170×L43×W24mm	Net weight: 130g	Dimensions: H180×L43×W24mm



Sonic Face Cleanser
• F62Y

S-Cleaner



S-Cleaner

Double-effective Brush Head
Double Cleanliness

01. Customized DuPont antibacterial brush head

0.05MM flexible DuPont brush head customized for providing a comfortable body feeling and completely pore cleansing. Antibacterial design, being capable of standing upright, natural air-drying, dirt resistance, being of more durable.



02. Silicone gel brush head with patented sonic motor

Silicone gel brush head with patented sonic motor is effective for cleansing without stimulating to skin. Edible hydrophilic silicone gel material, safe to health and environment-protected, soft feeling and effective cleansing, removal of dead skin. Division design of coarse and fine particles is considerate of the skin in different facial areas.



03. Much more effective on skin cleansing by use professional sonic motor

Magnetic levitation design for small internal friction, soft host vibration, long-lasting service life, easy for gripping by hand, subtle vibration to hand, direct output of brush head power and effective cleansing.



04. Intelligent partition reminder

S-Cleaner has a function of intelligent partition warning to remind the customer by a sound to avoid excessive cleansing on the skin.

05. Good brush head gives a deep cleansing to protect the skin.

- Standard cleansing mode
A common pattern for most people's daily cleansing.
- Soft skin mode
With silicone brush head for delicate, sensitive skin.



SPECIFICATIONS

Product Name: S-Cleaner	Built-in battery: 600mAh lithium battery	Material: ABS, PBT, silicone
Rated power: About 3W	Charging time: about 2.5 hours	Net weight: about 265g
Charging method: 5V=1A, Micro USB	Dimensions: H167.4×L75.6×W53.5mm	




Ion Face Cleanser
•F62



i-Cleaner
Ion adsorption
Double cleansing

01. Double-effect brush head, double cleansing, double efficacy


Only 1-2 minutes a day, easy to wash away make-up residue, hair particles, make the skin more bright, so that the pores are more transparent, far better than the hand-washing better cleansing effect, so that in the follow-up the skin care products are more easily absorbed.



i-Cleaner adopts patented sonic technology, with double-effect brush heads for different skins, and the skin surface is more single to clean; the ion derivation function is synchronized, the pores can be absorbed, the deep cleansing, gentle somatosensory does not hurt the skin, and the net shines.

02. Ion-Import function, one machine dual purpose

i-Cleaner has built in ion-import function for dual use. After cleaning the face, replace the dedicated introducer, sonic massage with iontophoresis, which can help the nutrient absorption of skin care products and avoid the problem of bacterial growth caused by surface residues.



03. Two brush heads, three modes, suitable for all skin types

- **Deep cleansing mode**
Suitable for non-sensitive skin with strong sebum secretion and regular makeup.
- **Standard cleansing mode**
A common pattern for most people's daily cleansing.
- **Soft skin mode**
With silicone brush head for delicate, sensitive skin.

SPECIFICATIONS

Product Name: i-Cleaner	Built-in battery: 600mAh lithium battery	Charging time: about 2.5 hours
Rated power: About 3 W	Dimensions: H167.4×L75.6×W53.5mm	
Charging method: 5V=1A, Magnetic suction	Material: ABS, PBT, silicone, Stainless Steel	

01. Are you facing these problems?




Black headed
Strawberry nose problem



Large pores
Skin is not delicate



Vigorous oil
Skin is a "big oil field"



Cosmetics residue's
hard to clean thoroughly.

02. You need a nanoTime Beauty ultrasonic exfoliation instrument!

This product uses ultra-high frequency 28000Hz ultrasonic technology, leading the market of general 24000Hz products, specializing in facial skin care and research and development, higher-end, more leading.

- ① 28KHz Professional Ultrasonic Technology
- ② Ultrasonic exfoliation
- ③ Ultrasonic cleaning pore
- ④ Ultrasonic + micro-current pulling and shaping
- ⑤ Ultrasonic importing assisted absorption
- ⑥ Ultrasonic energy rejuvenates skin



Ultrasonic
cleaning pore




Ultrasonic importing
assisted absorption




EMS micro-current
pulling and shaping

03. Deep cleansing by Multi-effect Integrated Ultrasonic

Ultrasonic thermal energy helps to open the pores. Ultrasonic function carries out 28000 times of "smashing" movement every second, shattering the stubborn blackhead of the pore and accompanied by the spray of toner water, ION exporting synchronization further absorbs pore dirt, triple function synchronization makes the skin more permeable.



04. Dual cleansing functions of ultrasonic and ION make skin transparent and white



05. LCD screen display is easier to use



✦ SPECIFICATIONS

Working voltage: 3.7V ---	Charging time: About 2 hours	Standby current: Less than 50 uA
Operation mode: W Peeling, W Moist, W Lift	Net weight: About 116g	Consumed power: about 2.5 W
Dimensions: L45 × W25 × H178mm	Charger requirement: : Output 5V ---1A	Material: ABS, Stainless Steel
Keys operation: Long press the power button switch machine, short press to switch mode		



Ultrasonic Exfoliation
Instrument

•F68



scriber

Ultrasonic cleansing keratin,
Micro-current pulling V-face

nanoMix-H

A real cosmetics sprayer
special for lotion
toner, essence



nanoMix-H
•F92

01. Atomized evenly, deep nutrition

nanoMix-H atomizes skin care products evenly, makes it penetrating pores rapidly to foundation of skin in most mist statue, high efficiently absorption.

Supplement
both water
and nutrients

Light and
portable

Usb charging
easy to carry

03. 13° spray angle, 300mm spray distance, spray more evenly

Strong atomization power, tens of thousands of tests verified that spray is more homogeneous and not wasteful (it is recommended to spray about 150mm above the face).

02. An era of atomization makes skin care easier than ever

Fast absorption with a gentle spray, product-saving, hydrating without messing up your makeup and the effect is several times better than hand applied. All of above advantages make skin care easier than ever.

Fast absorption
of nutrients

Save 50% amount
of skincare products

Hydrating without
ruining your makeup

Light and portable

Help with hydrating
and long-lasting makeup.

Usb charging
easy to carry

04. Atomization absorption, scientific skin care

- Pressurization technology has strong atomization force, which can atomize 2000mPa. 5 viscous water like skin care products.
- Apply to toner, essence and other conventional skin care products, scientific skin care effect is better.

purified
water

Toner

emollient
water

Skin
water

Essence

05. With polypeptide muscle repair solution, the effect is more outstanding

Nanomix-H portable skin invigorating instrument combined with NanoTime Beauty polypeptide muscle repairing nutrient solution can make skin absorb quickly and activate skin.



SPECIFICATIONS


Atomized amount: About 3ml/min	Battery Capacity: 7.4V, 450mAh	Dimensions: W42 × L47 × H161mm
Spray Distance: 200-300mm	Full charge time: About 2.5 hours	Net weight: About 188g
Tank volume: 6ml	Rated power: About 3.5w	Material: ABS, PC

01. A skincare artifact with new technology that can atomize all kinds of toners

Draw on the principle of medical atomizing treatment, it can evenly atomize various aquiform skin care products like toner, lotion, essential water, restraining water and smoothing toner to take moisturizing and nutrition.



02. A skincare artifact with new technology that can atomize all kinds of toners



Micron atomizing skin molecules penetrate pores rapidly until the end of muscle in aerosol mixed state, enhances skin absorbing capability.

03. Technology optimization, more stable and reliable

- ① More scientific water tank design
- ② More optimized gas channel
- ③ Higher adaptability of new materials



nanoMix-A
Portable Deluxe Skin
Activating Apparatus

SPECIFICATIONS

Atomized amount: About 3ml/min	Battery Capacity: 7.4V, 450mAh	Dimensions: L44×W49×H170mm
Spray Distance: 200-300mm	Full charge time: About 2.5 hours	Net weight: About 201g
Tank volume: 6ml	Rated power: About 3.5w	Material: ABS, PC

nanoMiX-E

Atomize essence
Better absorption



nanoMiX-E
· F92E



01. Better absorption through atomization

Supplement
both water
and nutrients.

Light and
portable


Use charging
easy to carry

02. Scientific skin care with atomizing absorption

Pressure boost technology can strengthen the atomizing function so that it could atomize aquiform skincare products of high viscosity (up to 3000mpa.s).

It applies to ordinary aquiform products such as toner, essence, sk2, etc, and provides a better skincare effect.


Micron-size droplets, Rapidly penetrate into the skin Easier to absorb, more effective skin care.



03. New design, easier to clean

If the misting of thick liquid fails or the mist volume is obviously small, the following steps can be referred to clean the nozzle.

- After removing the nozzle, press the convex pressure bar at the bottom of the nozzle with two fingers (as picture 1), and then press it down (as picture 2) to take out the nozzle accessories.
- Rotate the backup gear of the accessory anticlockwise and pull out the needle fitting slightly backward (as picture 3).
- Rotate the stainless steel part at the front of the fitting clockwise to remove the nozzle baffle (as picture 3).
- Wash all parts with clean water to clear the blockage. (please handle with care, do not clean and dredge with hard objects, so as to avoid damage to precision components).
- When installing the components, please strictly follow the above steps in reverse order.



SPECIFICATIONS

Atomized amount: About 3ml/min	Battery Capacity: 7.4V, 700mAh	Dimensions: L43×W50×H176mm
Spray Distance: 200-300mm	Full charge time: About 2.5 hours	Net weight: About 250g
Tank volume: 6ml	Rated power: About 3.5w	Material: ABS, PC



handy mist
sprayer

• F72

Handy 2s

Nano-atomization technology
for water replenishment
without dirtying makeup

01. Nano atomization technology, which makes the spray more delicate and achieves the purpose of replenishing water without dirtying makeup.



02. Intelligent design, slide cover spray, easy to use.



Fashion design
slide can spray



30 seconds timer, automatic
stop at that time



Small and portable
replenishing water at any time
and anywhere

03. Large capacity water tank

The capacity of the water tank is 8ml, which can be sprayed 8 times (30 seconds / times) with water.



04. Handy 2s Usage method



Step1. Open the water tank lids



Step2. Pick up the plugs of the water tank



Step3. Introduce clear water, cover the plugs and lids



Step4. Open the lid to spray

✦ SPECIFICATIONS

Product name: handy mist sprayer

Tank volume: 8ml

Spray time: one time 30 seconds

Battery capacity: 600mAh

Net weight: 77g

Dimensions: L48×W20×H115mm

Material: ABS, PC, Stainless Steel(Atomizer)



depiTime
Intensity hair removal

Ultra intense pulsed photons for
efficient and rapid hair removal

Super Photon
Dehairing Instrument

• B62E



01 Working principle of IPL hair removal instrument

This product uses a new generation of IGBT high-speed and high-energy photonic driver chip, then the IPL photonic energy is rapidly emitted through the epidermis and absorbed by hair follicle melanin, thus evaporating the nutrients of hair papilla, which can inhibit hair growth after several cycles.



Multi gear selections

Multi-levels can be selected for a variety of skin and all parts of hair removal



long life

500,000 extra-long service life, single/continuous hair mode is optional, enough for multiple people to use depilation



Intelligent protection

Intelligent launch protection to prevent misoperation, safety and durability



In the first month

Hair growth will slow down, hair removal is done once every week or two.



In the second month

depilation was performed every two weeks. Hair will gradually become thinner and thinner.



Starting from the third month

It can be operated once a month. Then it can be used selectively according to the growth state of hair.

Notice: Depending on the difference of human body constitution, the effect of hair removal will be different. After hair removal, it will not grow again in a few years. In special cases, it can be used repeatedly to gradually inhibit the growth of hair.



SPECIFICATIONS

Material: ABS, PC, Quartz (lamp tube), Silicone

Irradiation area: 30×10mm

Rated Power: 48W

Specification of adapter: AC 100V~240V 50/60Hz 12V=4A

light source: Xenon Flash Lamp

Net weight: About 300g

Preset Lifetime of Lamp Tube: About 500000 irradiations (at level 1) Dimensions: L72×W39.5×H204mm



firmAbs

Professional interference
wave products derived from Japanese
medical and nursing technology

MicroCurrent Interference
Wave Abdominal
Beauty Instrument


• B88



01. This product uses 4500Hz micro-current transported wave from Japan's professional medical and nursing field. It has no obvious pain, but also can pass through fat directly to the deep muscles. With 2-100Hz low-frequency interference current, it massages muscles vigorously to achieve high-quality EMS massage.



02. Dual independent EMS micro-current, Interaction of Four-Electrode Interference Wave Signals, Professional technology originated from medical and nursing instruments in Japan.



03. Professional interference wave micro-current has no stinging pain, which makes the body feel more comfortable; double waveform massage makes the function deeper; one minute massage is better than dozens of sit-ups, and can easily exercise, shape and beauty abdomen.



04. The main machine can be disassembled, so as to avoid the disadvantage that the product can not be used when the patch is out of order, which is convenient and sanitary, and the cost of long-term and multi-user is lower.



SPECIFICATIONS

Charging time: About 2 hours	Battery Specification: 3.7V, 180mAh	Charger requirement: 5V=1A
Dimensions: L55×H43×W16.5mm	Consumed power: About 1.5W	Standby current: ≤30uA
Material: ABS, Silicone, Copper (electrode)	Net weight: About 24g	



Tulipcurler

Simple and Fast Variable
Hair Style Anti-scald design,
hairdressing safer



Anti Scald Automatic
Hair Curler

• B53


01. Two-way automatic curling, simpler and faster to use

Two-way 360 degree rotation, easy automatic curling, you can make beautiful curls in both directions clockwise and counter-clockwise. Hair clips are soft, not winding, not pulling hair, not injuring hair, simple and convenient, which allows you to change your hairstyle in seconds.



02. Full anti-scald design, which allows beginners to use for the first time without fear of scald

Innovative tulip shape, unique anti-scalding structure design, even close to the face when curling hair will not occur scald, so that novices do not have to worry about the risk of scald when using ordinary curling bar products!



03. New PTC ceramic heating chip, more precise temperature control, more protection for hair

This product uses a new PTC ceramic heating chip, which is different from the common heating products on the market. PTC heating is faster and more sensitive. It can avoid the damage of long time high temperature or short high temperature to the hair, and the temperature is controlled in more precise. It can automatically stabilize when the set temperature is reached, and avoid the damage of hair caused by sudden high or low temperature.




04. It can nourish and repair hair synchronously with special hair care essential oil.

With the use of nanoTimeBeauty special high-temperature hair care essential oil, nutrients are more easily penetrated into hair scales open under the action of high temperature, to provide deep water and nutrients for hair, and to repair hair from the root.

05. Intelligent reminder function can avoid hair damage caused by too long perm time

This product will automatically emit buzzing tips when perming, to ensure that the perfect shape is made without harming the hair due to the long stay time, and to provide more protection for the hair of the changeable shape.



06. Three-levels temperature control, suitable for all kinds of hair

This product has built-in three-levels temperature control function to ensure rapid perm styling, while avoiding excessive temperature damage to hair.

220°C 209°C 180°C

SPECIFICATIONS

Rated voltage: 100-240v	Rated frequency: 50/60Hz	Rated power: 33-400W
Material of heating pipe: PTC	Material: ABS, Aluminum Alloy	Net weight: About 395g
Dimensions: L65.5×W65.5×H260 mm	Size of heat pipe: 28mm	working temperature: 180°C-220°C



An instrument for making
disinfectant water

•B91X



Rapid generation of disinfectant
water, simple convenient and
more affordable

01. Are you worried about the trouble caused by viruses and bacteria?

- Ordinary disinfectant must be diluted scientifically to use, alcohol is flammable and explosive, expensive and unsafe, more inconvenient?
- Household disinfectant water production instrument, Japan's electrolysis technology, three different concentrations, no need to dilute, ready to use, simple, convenient, safe and reliable.



02. Japan's new electrolysis technology, rapid generation of disinfection water, home disinfection good helper

Using the 4th generation electrolysis technology in Japan, sodium hypochlorite disinfection water can be produced quickly by salt and water.



Get disinfection
effect quickly



Quick access to
disinfection water



safe and 99.9%
sterilization-effect



Ultra low cost



03. Three concentration modes, without dilution, can be used with preparation, simple and convenient, No residue side effect, safe and reliable

- The built-in three electrolysis modes can quickly produce three kinds of disinfectants with different concentrations: low, medium and high. It can be used in a variety of home disinfection scenarios, such as tableware, kitchen, living room, bedroom, toilet, or baby products, fruits and vegetables. It is easy to use without dilution.
- CE RoHS FCC Certification, sodium hypochlorite disinfectant is safe, no irritation, no residue, no side effects, sterilization and disinfection at home, safe and reliable.



SPECIFICATIONS

Bottle Capacity: 270ml	Dimensions: L92×W72×H231mm	Concentration: L: 100-200mg/L
Consumed Power: 5W	Net Weight: About 135g	M: 200-300mg/L
Working Voltage: 5V=	Material: PP	H: ≥400mg/L



UVC destroys DNA to kill bacteria, and photocatalyst technology realizes deodorization and antifouling



Portable UV
Air Purifier

•H04

01. Coronavirus is rampant, working from home is a lot of trouble?

How to effectively kill bacteria and viruses at home / office?
Common disinfectant is expensive, unsafe and inconvenient?



02. Photocatalyst sterilization synchronization, double effect disinfection more thorough

Photocatalyst technology can produce positive oxygen in the inner chamber under the action of ultraviolet rays, which can simultaneously kill bacteria and viruses. It has a dual effect, not only disinfection and sterilization, but also remove odors.



03. 253nm specific UVC sterilization and disinfection, the official recommended disinfection tools

The Portable UV Air Purifier adopts 253nm ultraviolet disinfection technology with specific wavelength which can kill all kinds of viruses and pathogenic bacteria.



04. Full airtight design, ultraviolet radiation does not leak, people can be more safe on site

No irritation, no residue, no side effect, certified by CE FCC RoHS ETL, efficient sterilization, bedroom, kitchen, toilet, baby room, pet shed, desk Safety disinfection is more assured.



05. Vortex inhalation, automatic air circulation, 24-hour continuous disinfection is more convenient



06. Portable design, plug-in can be used at any time, keep every room healthy and fresh



◆ SPECIFICATIONS

Product name: Portable UV Air Purifier	working voltage: 100-240V	Working noise: ≤60dB
Dimensions: L65×W64×H186mm	Rated power: 4.5W	UV light source: 3W/253.7nm



Aqua Perfect

Anti-oxidation
Beautifying And Anti-aging



AquaPerfect
•B91

01. High hydrogen content

Under room temperature, the theoretical concentration value of the soluble water element in water may be up to 1000ppb. If using simple tester, there is great error. It is not possible to be normal 1000 ppb under the room temperature.




02. Special degerming Mode

The Aqua Perfect Hydrogen Water Anti-oxidation Instrument adopts unique Japan's Hirose patented sterilization mode in this industry to resolve cleaning problems, more hygienic and well.



03 Supporting the addition of hydrogen water purifying agent

The hydrogen water produced by Aqua Perfect Anti-oxidation Instrument may relieve accumulation of the lactic acid and remit fatigue and muscle soreness arising from the exercises, relax your muscle and skin, rapid water replenishing while scavenging free radical.



04. Beautifying and anti-senility

If the hydrogen water produced by the Aqua Perfect is taken orally, it may rapidly enter into the human body's blood circulation system to eliminate free radicals and improve skin conditions from inside to outside; the hydrogen water produced by the Aqua Perfect has anti-oxidation, anti-inflammatory and anti-anaphylaxis effect, effectively relieves the aging effects of the free radical to skin to improve dry skin, wrinkle and discoloration problems.



05. Improvement of human body's metabolism function, bring health and beauty

Due to the molecular properties of water with hydrogen element produced by Aqua Perfect, the grease may be dissolved. Therefore, it has good health care function to three high problems.



water



Hydrogen

The test proves that the water with high-concentration hydrogen water may dissolve grease; in the ordinary water, the grease is floated on the top of water.

✕ SPECIFICATIONS

Product name: AquaPerfect

Product Model: B91

Volume: 300mL

Water temperature range: 5-30°C

Battery: 1200 mAh

Adapter: 5V/1A

Hydrogen concentration: can reach 600 ppb in 5mins

Net weight: About 226g

Dimensions: #9 Φ72*H198mm

Material: ABS, PC, Silicone, Titanium alloy (electrode)

01. The specially customized photon cleaning instrument is matched with the recommended peroxide bright white gel by the dentist. The blue light with certain wavelengths is used to catalyze the bright white active factor.



02. Whitening gel pens for special quantitative filling are simpler and more convenient to use.



03. One-to-three USB interface, which can be directly powered by mobile phones, is more convenient to use.



Photon Dental
Cleaner

·IT01

04. The special bright white gel contains safe and reliable teeth cleaning ingredients. The active whitening factor helps to remove stains on the teeth surface, while protecting your enamel.




Active
ingredient



Specific
wavelength



Bright
white gel



Photon Enhanced Active
Factor Action Dual effects
make teeth brighter and whiter

◆ SPECIFICATIONS

Product Name: Photon Dental Cleaner	Power supply interface: Micro USB	Rated voltage: 5V=
Consumed power: About 1.2 W	Current: About 240uA	Net weight: About 26g
Dimensions: H65×L37×W24mm		



Magnetic deep assisted absorption+ red light to skin rejuvenation Easily eliminate eye skin troubles



Eye beauty
•F51

01. Do you face eye skin problems?

With the growth of age, the problems of eye skin have not been solved, even more and more serious. Black circles, fine lines, pouch...A variety of eye problems continue to emerge.

Is it difficult to absorb nutrients in the eye skin? No matter how expensive eye creams and eye essences are used, not only is it useless, but even fat particles have grown due to excess nutrition?



02. Three functions to easily solve eye troubles



- **Magnetic induction**
in multiple layers, and the magnetic energy accelerates the introduction, so that the skin around the eyes can absorb nutrients more efficiently.
- **Sonic**
Sonic massage vibration relieves eye fatigue, promotes eye circulation, and reduces eye line lines.
- **Red Light**
Red light skin rejuvenation, restore the under-eye skin, dilute dark circles, and improve eye problems.

03. Well-designed product details

Magnetic massage head

Integrate LED phototherapy, sonic massage, magnetic induction function electrode in one, precise technology, just to make you more beautiful.



Smart touch design

No key operation, hold the electrode workinghead of the machine body and contact the skin to start working.



Ergonomic structure

Scientific inclined plane design for facial, eye and jaw massage. Light weight design, easy to carry, no space occupied.



✦ SPECIFICATIONS

Product name: Eye beauty

Product model: F51

Battery Specification: AA battery

Consumed power: 0.3W

Net Weight: About 28g



Visible clean skin,
precise sniping
of "blackheads"



visual blackhead
instrument

•F31



01. Three functions assist in skin care, making pore cleansing easier

Viewing mode

HD camera, with LED fill light, 30 times magnification display, so that blackhead have nowhere to hide.

Blackhead mode

three levels of suction power can be adjusted, and different suction heads can be selected according to different skins, which is perfect to suck blackheads without hurting the skin; the entire operation process is visible, more convenient and cleaner.

Blue light mode

415nm wavelength care blue light, which acts on the skin surface, induces the apoptosis of Propionibacterium, plays a bactericidal effect, effectively inhibits the secretion of sebaceous gland oil, reduces pore clogging, and effectively inhibits the production of blackheads and acne.

02. Three suction heads for different parts

Oval head for nose



The big round head fits the nose



Small round head suitable for nasolabial groove area, or stubborn blackhead nursing



03. 3D curved nozzle design

Fits the curved surface design of the face, the outer edge is round and smooth. It can better fit the skin and strengthen the suction, but it is not easy to scratch the skin. It is soft and smooth and sucks out the dirt.



04. Baby pacifier grade silicone

Flexible, not easy to deform, gentle to the skin, reject the damage caused by the plastic hard nozzle.



05. Smart mobile reminder to avoid skin damage

In the blackhead mode, the blackhead app can intelligently identify whether excessive suction occurs in an area, and prompt you to move the blackhead instrument through the screen in time to avoid damage to the skin.





Moisturize your skin at any time, Enjoy technology skin care



Hydrogen Water Skincare Atomizer
•F97



01. What is hydrogen water ?

Hydrogen water is saturated hydrogen water produced by electrolysis of high quality electrode. Its release of micron molecules has strong diffusivity, which can deeply moisturize skin, anti-aging.

02. Moisturize the skin anytime and anywhere, and make the skin relaxed all day

Moisturizing and beautifying

Through the built-in titanium electrode sheet, electrolyze water rapidly to produce ultra-fine hydrogen ions, penetrate into the cells, eliminate the active oxygen in the body, help the skin metabolism, and easily achieve the effect of anti-aging and beauty.



Superfine mist

Hydrogen molecules diameter of water element is 10 to 150 000 of human cells. It can directly enter human cells and deliver nutrients to the bottom of muscles. This can effectively moisturize the skin and fade the fine lines to plump and beautify the skin.

Smart and portable

Compact design, one-handed control, can be carried when traveling. Moisturize the skin anytime and anywhere.



SPECIFICATIONS

Product Name: Hydrogen Water Skincare Atomizer	Product Model: F97	Tank Capacity: 17ml
Spray Volume: 20-40s/3mL	Rated Power: 4.6W	Working Time: ≥40 mins
Battery Capacity: 1200mAh	Material: ABS	



Precise Effective
Skin Beauty



Freezing Point
Photon Hair
Removal Device
•B62C

01. NanoTime Beauty Freezing Point Photon Hair Removal Device, give you a new experience of cold and comfortable hair removal



02. 6 major advantages, comfortable and safe to hair removal

1. Sapphire automatic ice compress nursing, ice compress and depilation synchronization



2. High energy color light to soften skin, depilation more tender skin



3. 500000 extra-long service life, shared by the whole family



4. High energy continuous hair mode, complete the whole body hair removal in 10 minutes



5. 6 gears can be selected to meet different hair removal requirements



6. Intelligent protection design to prevent misoperation



03. Sapphire + semiconductor refrigeration technology = 10°C ice sensation throughout

Powerful freezing point hair removal, no pain and no damage to the skin, easy to smooth skin



04. Double mode to choose as you want, making hair removal more convenient

Manual mode small area hair removal flexible and free
Continuous flash mode large area hair removal effortlessly





Straight hair
comb + curling bar,
one machine for dual use



Intelligent multifunctional
hairdressing instrument

• B56

01. Professional styling hairdressing instrument, easy to make a variety of styles

- Dual use of rolling and straight
- Available in different curls
- Full anti-scalding design
- High-definition display
- Intelligent temperature control
- Ion hair care



02. Three curly hair diameters + straight hair comb free choice

19mm diameter, lazy small wave, full score of vitality. 25mm diameter, fashion medium volume, create a beautiful temperament. 32mm diameter, the charm of big waves, shining confidence in the crowd. Straight hair comb, easy to smooth the hair quality changes visible to the naked eye



03. 360° left and right rotation control axis

- For personal use ✓
- For other use ✓
- Involution ✓
- Eversion ✓

04. Anti-scald straight hair comb

The comb teeth adopts heat insulation and anti scald design, easy to comb hair, easy to use, and more rapid.



05. Four levels of precise temperature control, suitable for all hair types

- 220°C
- 200°C
- 180°C
- 160°C

06. HD display

Real-time monitoring of heating temperature and heating time.

07. Enjoy "Prince with long hair"

Extended heater Refuse to roll half of the embarrassment, and waist long hair can also be in place.



New artifact of disinfection and epidemic prevention, black technology of fresh air



Air Disinfector
·H05

01. Triple cycle of disinfection and sterilization, which has a new definition of air purification

First Physical Filtering

H05 air disinfector adopts a professional HEPA filter screen with Level F8 filtration rate, which can conduct physical filtration for various PM2.5 particles in the air, more reliable and safer.



Third UVC inactivation, photocatalyst synchronization

H05 air disinfector adopts medical grade 253nm + 185nm dual band customized UVC light source and photocatalyst sterilization technology from professional medical scene and laboratory level, which can kill all kinds of bacteria and viruses, and the killing rate can reach 99.9% according to the report of authoritative testing agency. At the same time, its photodegradation combined with photocatalyst will produce trace ozone to further catalyze the decomposition of formaldehyde, phenylammonia, peculiar smell, etc.

02. Full protection design, more safe disinfection and sterilization




Second Electrostatic adsorption, ozone killing Bacterial

H05 air disinfector adopts laboratory level electrostatic adsorption technology to further absorb various dust impurities, bacterial viruses, etc. in the air on the basis of physical filtration. Meanwhile, it produces a certain concentration of ozone to kill bacterial viruses and decompose harmful substances such as formaldehyde, benzene, ammonia, etc.



03. Replaceable light source, replaceable strainer longer service life



04. Eddy current active circulation sterilization, no dead angle, more sustainable



Disinfection and
sterilization easily
protect family health

Air Disinfectant
•H08



01. Quadruple cycle disinfection and sterilization, a new definition of air purification, multiple protections for family health

HEPA Filtering PM2.5+Electrostatic Adsorption / Ozone Disinfection+ UVC Photocatalyst Photolysis Disinfection+ Disinfection Water Washing. It can solve the problems of pollen allergy, bacteria breeding in kitchen / toilet, difficult to remove odor and formaldehyde carcinogenesis



02. How is this product sterilized?

The First. Physical Filtering

It's adopts a professional 8-layer HEPA filter screen, which can conduct physical filtration for various PM2.5 particles in the air, which is safer and more reliable.




The Third. UVC inactivation, photocatalyst synchronization

H08 air disinfectant adopts medical grade 253nm + 185nm dual band customized UVC light source and photocatalyst sterilization technology from professional medical scene and laboratory level, which can kill all kinds of bacteria and viruses. And catalytically decompose formaldehyde, benzene ammonia, peculiar smell, etc.

The Second. Electrostatic adsorption, ozone killing

H08 air disinfectant adopts laboratory level electrostatic adsorption technology to further absorb various dust impurities, bacterial viruses, etc. in the air on the basis of physical filtration. Meanwhile, it produces a certain concentration of ozone to kill bacterial viruses and decompose harmful substances such as formaldehyde, benzene, ammonia, etc.



The Fourth. Disinfecting water flushes the air

It uses electrolyzed salt water to produce sodium hypochlorite solution to wash the air again, kill the possible residual bacteria and viruses, and further oxidize and decompose formaldehyde, aniline and peculiar smell. When the air is washed by hypochlorite solution, the residual ozone molecules will be dissolved in water and decomposed into oxygen when meeting water, so as to completely avoid the influence of ozone escaping on respiratory tract.



03. Smart sensor function

Automatically recognize the air quality, and intelligently start the corresponding aldehyde removal mode, deodorization mode, disinfection mode, and standard mode.