







ANHUI JUHUI TRADING CO.,LTD

- admin@ahjhmask.com
- 86-556-5800966
- Room 313-2,3rd Floor ,Bldg A1 ,Zhumeng New Zone,No.188,Wenyuan Road,Anqing,Anhui,China

**Professional and Customer-Oriented** 

安徽聚辉贸易有限公司 ANHUI JUHUI TRADING CO.,LTD

# Company Introduction

With over 15 years of experience Anhui Juhui Trading Co., Itd is your reliable, expert and flexible partner in disposable, protection and hygiene products for the medical, industrial and cosmetic sector.

You and your customers are at the heart of our business

We place great emphasis on intensive interaction with our domestic and overseas customers.

We are motivated by your needs and wishes; the benefit to you is informative and solution-oriented support.



# Juhui—Expert, Professional and Customer-Oriented

- Reliable service
- Solution-oriented & fair
- Detailed consultation
- Customized product design,
- packaging and storage extensive product range
- Tried-and-Tested quality
- Modern warehouse logistics



# Juhui – Flexible, Efficient & High-performance

- High supply capability
- Fast delivery
- Optimal price/performance ratio
- Quick response
- Friendly & experienced team

#### Juhui-Superior Quality, Extensive Mask Range

- Disposable non-woven face mask
- Folding shape protective respirator face mask
- Butterfly shape protective respirator face mask
- Cup shape protective respirator face mask
- Fish shape protective respirator face mask
- Duckbill protective respirator face mask
- Printing protective respirator face mask
- Anti dust/poisonous protective respirator face mask

# Product catalogue

03-04	Disposable Non-woven Face Mask
05-06	Folding Face Protective Mask
07	Butterfly Shape Protective Face Mask
08	Fish Shape Protective Face Mask
09	Duckbill Shape Protective Face Mask
10	Gauze Protective Face Mask
11-12	Cup Shape Protective Face Mask
13-14	Printing Protective Face Mask
15-16	Anti Dust/Poisonous Protective Face Mask
17-18	Dust Sponge Protective Face Mask

#### Certificate









FFP2

**NIOSH** 

CE

SGS

#### Our Quality Management System ISO 9001:2015 Certified

- The continuous improvement process is an essential part of or internal processes.
- Our explicit goal is being efficient and effective in order to offer you high quality at fair price, in accordance with your needs.
- We pay high attention to the quality of our products.
- In this context, careful supplier selection is essential. In addition to the relevant certifications, we base our decisions primarily on reliability, quality, and innovative products, as well as the observance of internationally recognized conventions governing working conditions and standards.





Size :17.5cm\*9.5cm

Color :Pink,white,blue,green,

Black, Gray

Style :Tie on and ear loop

(round or flat elastic)

Packing: 50pcs/box, 40boxes/carton or

as customers requirement.

Feature: High fluid respiratory protection,

easy breath, anti-flu, high filtration

Usage: Hospital, foodprocessing,


beauty salon, electronic industry, etc.







# Disposable Non-woven Face Mask





• Fashion appearance

Professional design

High-quality materials

Import Japanese Packaging

Ordinary English Packaging



 Four Layers of Activated Carbon Mask



Korea Individual Packaging





ITEM NO	WEIGHT	LAYER	FILTER
JH7301	18+18+25g	3Ply	BFE>99%
JH7302	20+20+25g	3Ply	BFE>99%
JH7303	30+20+25g	3Ply	BFE>99%
JH7304	20+20+25g	3Ply	BFE>99%
JH7305	20+35+20+25g	4Ply	BFE>99%





Layer :5ply,6ply

Color :Pink,white,blue,black,gray.
Style :Latex-free elastic headband

and ear loop

Packing: 25pcs/box,500pcs/ctnor30pcs/box, 300pcs/ctnor40pcs/box,800pcs/ctn.

Feature: Anti-dust , Anti-fog, Anti-haze,
Anti-PM2.5, Anti-germs.

Usage: Home-care, hospital, clinic, outdoors.







# Folding Face Protective Mask

## **Protection Efficiency**

1st layer	* * *	PM 10 PM 3 PM 2.5 PM 0.3
2nd layer		PM 3 PM 2.5 PM 0.3
3rd layer		PM 2.5 PM 0.3
4th layer		PM 0.3



# Packaging









ITEM NO	Belt Type	Packing	Filter
JH8501	Ear-loop/Headband	10/25/60pcs/box	KN95
JH8501V	Ear-loop	25/30/40pcs/box	KN95
JH8502	Ear-loop/Headband	25/50/60pcs/box	KN95
JH8502V	Ear-loop	25/30/40pcs/box	KN95
JH8602	Headband	25/50/60pcs/box	KN95
JH8602V	Headband	25/30/40pcs/box	KN95





:5ply,6ply Layer

:Pink,white,blue,black,gray. Packing: 25pcs/box,500pcs/ctn or 30pcs/box,300pcs/ctn or 40pcs/box,800pcs/ctn.

Style :Earloop

Feature: Anti-dust, Anti-fog, Anti-haze,

Anti- PM2.5, Anti-germs.

Usage: Home-care, hospital, clinic,

outdoors.



:4ply,5ply Layer

(available with active carbon)

Color :White,blue,black,brown,pink Packing:2pcs/bag,30pcs/box,600pcs/ctn;

2pcs/bag,20pcs/box,400pcs/ctn

Style :Latex-free elastic ear loops

Feature: For protection against dust, germs, PM2.5, fog and haze.

Usage :Suitable for persons

wearing glasses















**Butterfly Shape Protective Face Mask** 

#### Fish Shape Protective Face Mask

#### **Material**













- Non woven out layer, fluid resistant and flame retardant.
- High efficiency meltblown filtrtion layer, At least 95% filtration efficiency.
- 3 Needle punched cotton layer
- Thin inner layer, soft and skin-friendly

#### **Material**

- 11st layer: non-woven polypropylene ≈ 48 gsm
- 2 2nd layer: hot-blast nonwoven ≈ 35gsm
- 33rd layer: melt-blown nonwoven ≈ 30gsm
- 44th layer: non-woven polypropylene ≈ 48 gsm









Elastic straps: spandex

# **Specification**

ITEM NO	Valve	Protective Level	Color
JH8501-2	Without	N90/N95/N99	Black,Gray,White,Pink,Blue
JH8501-2V	With	N90/N95/N99	Black,Gray,White,Pink,Blue

ITEM NO	Valve	Protective Level	Belt Type
JH8401	Without	KF80/KF94	Earloop
JH8401V	With	KF80/KF94	Earloop







Color :White,blue,grey,black

Packing: 300pcs/ctn or 600pcs/ctn

Style :Headband

Material: Non-woven fabric, heat-sealing

Sponge, Melt-blown fabric.

Usage : Appropriate for situations where exposure to blood and/or bodily

fluids is not a risk.





:9ply,12ply,16ply,24ply

Color :White

Packing:500pcs/ctn or 800pcs/ctn

or 1000pcs/ctn

:Tie on and ear loop

Material: Non-woven fabric/gauze Feature: Reusable and washable. Usage :Hospital,food processing,

electronic industry, etc.





#### **Duckbill Shape Protective Face Mask**







#### Gauze Protective Face Mask










How to use





Divide up and down Make mouth of mask downward Wear the mask of the

of the mask by hand and the headband overhanging face and pull the headband





should be on the nack





The bottom band of mask Top band of the mask should Press the nose dip and be on the top of head

**Material** 



# Specification

ITEM NO	Valve	Protective Level	Color
JH8301	Without	N90	White, Blue, Gray, Black
JH8301V	With	N95	White, Blue, Gray, Black

ITEM NO	JH7009	JH7012	JH7016	JH7024	JH7501
Layer	9Ply	12Ply	16Ply	24Ply	Mechanism of masks
Packing	1000pcs/ctn	800pcs/ctn	500pcs/ctn	500pcs/ctn	1000pcs/ctn





:3ply,4ply Layer

:White,black,gray. Color

:Latex-free elastic ear loops Packing:10pcs/box, 200pcs/ctn or

20pcs/box, 400pcs/ctn;

Feature: For protection against dust,

germs, PM2.5, odor, fog and haze

Usage: Suitable for persons wearing glasses







Packaing

# **Cup Shape Protective Face Mask**

# Material











- Non-woven, fluid resistant and flame retardant.
- 2 MeltBlown fabric, at least 95% filtration efficiency.
- 3 Acupuncture cotton, soft and skin-friendly.



## Specification

ITEM NO	Valve	Protective Level	Color
JH9301	Without	N90/N95	Black,Gray,White
JH9301V	With	N90/N95	Black,Gray,White
JH9101	Without	N90/N95	Black,Gray,White

#### Application













Material: Non-woven fabric, heat-sealing

sponge,melt-blown fabric.

Layer: 5ply,6ply

Color : Animal, star, snow or customize

printing

Style: Ear loop
Packing: 40pcs/box,800pcs/ctn or

25pcs/box,500pcs/ctn or 30pcs/box,300pcs/ctn









Usage:Hospital,food processing,electronic industry,etc.

#### How to use



Face the mask side without the nose strip, make the nose strip is above the mask.



Put the mask on the chin.



Put the ear loop around your ears.




Put the finger in the middle of nose strip,
Pressing the nose strip from the middle to
the two sides, until it is completely match the
shape of nose.



Check the adhesion of masks before entering the

#### **Printing Protective Face Mask**



ITEM NO	Valve	Protective Level	Belt Type
JH8201	Without	N90/N95	Earloop
JH8201V	With	N90/N95	Earloop





Type: Head loop Color: Gray Material: Rubber, TPR Plastic Filtration rate: 95%

#### **Application**

Grinding, sanding, sweeping, sawing, bagging, or processing minerals, coal, Iron-ore, flour, metal, wood, pollen, and certain other substances.liquid or non-oil based particles from sprays that do not also emit oil aerosols or vapors.



#### **Feature**

- High quality rubber cover, good oil resistance, comfortable durable, suitable for long wear, easy to clean and maintain.
- Headband head shape design, more fit, strong and durable.
- Recyclable, replace the dust filter core. For protection against dust, germs, PM2.5, fog and haze.

#### Material



#### Anti Dust/Poisonous Protective Face Mask



Item No	Packing	Protective Level
JH670	1pcs/box,40box/ctn	N95
JH672	1pcs/box,40box/ctn	N95
JH674	1pcs/box,40box/ctn	N95





Material: Sponge, polyester cotton, other

materials are available

Type :Ear-wearing

Thickness:2.0mm,2.3mm,2.5mm

Packing:1pcs/vacuumbagor3pcs/vacuumbagor1000pcs/ctn,5000pcs/ctnor

bag or rooopes/ein,sooopes

customer's size.

Function: Anti-dust, Anti-fog, Anti-haze,

Anti-PM2.5, Anti-germs.

Feature:Breathable, comfortable,

skin friendly













## Wearing Steps







Wear the mask



Adjust the mask to fit your face

#### Notes

- Not applicable to hazardous or gas leakage
- If you feel skin uncomfortable, please stop using it immediately
- If your feeling of nausea or odour is umbearable, please stop using it immediately
- Avoid placing it in high temperature, high humidity or direct
- Wearing glasses may cause fog, please pay attention to such situation when driving

#### **Dust Sponge Protective Face Mask**



#### **Specification**

Item No

JH1001,JH1002,JH1003,JH1004,JH1005

Color

Red, Pink, Purple, Green, Blue, Yellow, White, Black, Gray.