



Hot/Cold Therapy
More Healthy Life



江苏英科医疗制品有限公司
Jiangsu Intco Medical Products Co., Ltd.
Add/江苏省镇江市大港新区烟墩山路77号
No.77 Yandunshan Road,Zhenjiang 212132
Tel/86-0511-83174088 Fax/86-0511-83174088
Http://www.intcomedical.com



Instant Cold Pack & Instant Hot Pack

Instant Cold Pack.....01
Instant Hot Pack.....02

Reusable Hot & Cold Products

Gel Pack.....04
Gel Beads.....07
Clay Pack.....08
Flexiable pack.....09
Massage Ball.....10

SPA Series

Moist Heat Wrap.....11
Magic Warmer.....12
Eye Mask.....13
Face Mask.....14

Disposable Warm Patch

Normal Heat16
Pain relief.....17

Home Care Series

Ice Bag.....19
Fever Cooling Patch/Reusable Cooling Pack.....20
Cooling Gel Mat/PCM Cooling Mat.....21
Cooling Gel Cushion.....23
Cooling Gel Pillow.....24
Cooling Scarf/Cooling Belt.....25
PVA Towel.....26

Mom & Baby Series

Breast Gel Pack.....28
Perineal Instant Cold Pack.....29
Infant Heel Warmers.....30


Promotion & Gift Series

Warm Toy/Magic Warmer.....32

『英科目标』 让英科成为
中国最优秀民营企业之一

『英科精神』 诚信正直 勤奋刻苦
专业精神 团队合作

Instant Cold Pack



First Aid/Sports Support/Home Care.
External application only.
For relief of pain and swelling caused by sprains, strains, contusions, minor burns, toothaches and insect bites.

Item Number	ICE 0506B	ICE 0609
Description	12.7*15.3CM	15.3*22.9CM
Weight	109G	237G
Material	PE&NYLON	PE&NYLON

Characteristic:


Single use cold application without the drippy mess of ice cubes. Needs no refrigeration, convenient and portable. Primary usage for outdoor activity injuries and first aid.

Usage:

Locate inner pouch. Squeeze firmly to rupture inner fluid pouch. Shake briefly to mix contents. Wrap in soft cloth and apply. Best if activated between 32 and 80 degrees Fahrenheit.



Instant Hot Pack



Instant relief for knocks, bruises, sprains and sports injuries. Relieve frostbite in cold season, or arthritis in wet weather. Improve blood circulation.

Item Number	HOT 0506	ICE 0609
Description	12.7*15.3CM	15.3*22.9CM
Weight	104G	247G
Material	PE&NYLON	PE&NYLON

Characteristic:

No heating required. Ideal for emergencies, work instantly. Convenient and portable.

Usage:

Remove the inner bag by hand, do not cut with scissor or a knife. Crush the inner bag for a few seconds to activate the contents. Apply the heat pack to clothing and not the skin.



Reusable HOT& COLD Products

Used either for cold or hot
Specially designed for parts of body
Use as flexible ice or soothing heat
Cold helps reduce pain and swelling
Heat helps muscle aches and stiffness

Reusable HOT& COLD Products Gel Pack

Cold therapy application

Reduce the swelling and pain caused by inflammation and relieve the discomforts associated with:
Joint & muscle injuries / Back & neck pain
Arthritis / Sinus & stress headaches
Post-surgical pain / Sunburn & other burns

Heat therapy application

Relax tight muscles and restore flexibility in a wide variety of ailments such as :
Muscle pain & soreness
Muscle spasms & cramps
Arthritis / Menstrual cramps
Chronic neck or back stiffness



			
Item Number	HC-GEL0507	HC-GEL0507	HC-GEL0510
Description	12*17.8CM	17*25CM	12.7*25.4CM
Weight	144G	375G	325G
Material	PE&NYLON	PVC	Fabric

			
Item Number	HC-GEL 0511	HC-GEL 0410	HC-GEL0812
Description	12*29CM	10.2*24.5CM	20*30CM
Weight	350G	175G	700G
Material	PE&NYLON	PVC compound knitted fabric	PVC compound spunlace

an INTElligent COmpany to partner with

www.intcomedical.com



Usage:

To heat, microwave or boil the pack by instructions.
To obtain cold therapy, freeze the pack for at least one hour.
Place the pack in the bag of the support for usage.
Individuals with problems in sensing heat or cold should use extreme caution with this product.
For external use only.

Characteristic:

Remain soft and pliable after freezing.
Different sizes and colors are available.
Reusable with elastic wrap.

Item Number	WRA 0424	WRA 0709A	WRA 0609A
Description	9*61CM	18*21.8CM	16*21.8CM
Weight	535G	300G	300G
Material	PE&NYLON	PE&NYLON	PE&NYLON
Wrap	Fleece\Neoprene fabrication	Fleece\Neoprene fabrication	Fleece\Neoprene fabrication

Item Number	WRA 0816	WRA 0609B	WRA 06010
Description	40*20CM	22*16CM	16.5*25.9CM
Weight	850G	320G	320G
Material	PVC	PVC	PVC
Wrap	Fleece\Neoprene fabrication	Fleece\Neoprene fabrication	Fleece\Neoprene fabrication



Item Number: WRA 0822
Description: 19.7*56CM
Weight: 850G
Material: TPU



Item Number	WRA-R2540	WRA-R4060	WRA-R5060	WRA-R6060	WRA-R7080
Description	R2540	R4060	R5060	R6060	R7080
Weight	92G	565G	795G	910G	1360G
Material	TPU	TPU	TPU	TPU	TPU

Item Number	WRA 0721	WRA 1014	WRA1319
Description	17*53.4CM	25.4*34.5CM	32*47CM
Weight	910G	1360G	2270G
Material	TPU	TPU	TPU

Reusable HOT& COLD Products Gel Beads



This kind of hot and cold packs have good performance therapy for pain, tension, body aches, swelling and bruises with nice appearance.

Innovative soft gel beads provide long lasting cold relief and mold to your body.

Colorful gel-filled beads and form-fitting shape provide soothing relief.

Instruction:

Hot Therapy—Heat in microwave for specified time
Heat for additional 5 seconds as needed
until desired temperature is reached.
Cold Therapy—Freeze for at least 2 hours.

Item Number	HC-GB7523	HC-GB1317	HC-GB2038
Description	7.5*23CM	13*17.5CM	20.5*38.5CM
Weight	96G	245G	425G
Material	PVC	PVC	PVC

Item Number	HC-GB2894	HC-GB4124	HC-GB2819
Description	28*9.4CM	41*24CM	28*19CM
Weight	135G	450G	350G
Material	PVC	PVC	PVC

Reusable HOT& COLD Products Clay Pack



It utilizes a synchronized blend of all-natural, Ceramic clay and non-toxic oils to create the ideal thermal transfer source.

It capitalizes on the efficiency of radiant energy to create deep penetrating, therapeutic pain relief.

It is safe, easy to use, and highly efficient in the treatment of a wide variety of pain discomforts such as arthritis, sports injuries, muscle pain, joint stiffness, menstrual cramps and headaches.

It's hot performance is at least double effective than regular gels.



Item Number	WRA-CLAY 0612	WRA-CLAY 0916
Description	15*30CM	24*41CM
Weight	670G	1400G
Material	PE&NYLON	PE&NYLON
Wrap	Fleece\Neoprene fabrication	Fleece\Neoprene fabrication

Reusable HOT& COLD Products Flexiable pack



For instant relief from aches and pains. To provide soothing warmth and comfort, place this product in a microwave and heat in accordance with the heating guide below. Please always ensure that the microwave is clean and the turntable can rotate freely.




Great for headaches, migraines, sprains and bruising. Simply place in the freezer and chill for at least three hours. This product remains flexible when frozen, making it easy to apply to any area of the body.

		
Item Number	FPP1260	FPS1260
Description	12*60CM	12*60CM
Weight	705G	705G
Material	Kint fabric	Kint fabric

Reusable HOT& COLD Products Massage Ball



Massaging your body can help with slight aches, pains and tension in your back, arms and legs. The aim of massaging is to relax and stretch the muscles, thus promoting circulation and a free flow of oxygen, which ultimately benefits your sense of well-being.

			
Item Number	MB-6	MB-6	MB-8
Description	Diameter 6CM	Diameter 6CM	17cm, Diameter 8CM
Weight	105G	105G	220G
Material	PVC	PVC	PVC



SPA Series Moist Heat Wrap



Moist Heat has been used for hundreds of years by many cultures in the treatment of arthritis, rheumatism, joint stiffness and muscle spasms.

Moist Heat is preferred because of ability to penetrate deeply and more quickly at lower temperatures than dry heat.

Instruction:

For hot use, put the pack flat in microwave oven to heat, the heating time is depended on the different sizes.

For cold use, freeze the pack at least 2 hours in freezer before use.



Item Number	MH 250	MH 850A	MH 850B
Description	10.3*58.9CM	18.5*99CM	17.6*75.2CM
Weight	250G	850G	850G
Material	Fleece /on your need	Fleece /on your need	Fleece /on your need

Item Number	GH 310	GH 570	GH 680
Description	11*56.5CM	12*55CM	17.8*53.3CM
Weight	310G	570G	680G
Material	Fleece /on your need	Fleece /on your need	Fleece /on your need

SPA Series Magic Warmer



Moist Heat has been used for hundreds of years by many cultures in the treatment of arthritis, rheumatism, joint stiffness and muscle spasms.

Moist Heat is preferred because of ability to penetrate deeply and more quickly at lower temperatures than dry heat.

Instruction:

Fold the metal disc inside the pack, it will start crystallization. When the crystallization is completed, it's ready to use.

To recharge, please place heat pack in boiling water with cloth. Heat it about 10 minutes or until the crystal is dissolved. It's reusable.



Item Number	MAG 0507	MAG 0711	MAG 0816
Description	13*19CM	19.2*28.8CM	40.8*21.5CM
Weight	232G	700G	840G
Material	PVC	PVC	PVC

Item Number	MAG 0304	MAG 0410
Description	9.5*6.5CM	10.5*26CM
Weight	60G	150G
Material	PVC	PVC

SPA Series Eye Mask

Eye Mask

Subside black eyes and wrinkles of eyes to keep eyes young and beautiful. Relieve and recover readily. Reduce the overburden of the brain, such as migraine, insomnia etc.

Instruction:

Place eye mask into the refrigerator 1 hour before use. Wear it about 15 minutes.



"Age can always be a Secret"



Item Number	EYE-GB 0307	EYE-GEL 0310	EYE-GB 0310	EYE-GB 0309
Description	18.9*6.8CM	24.5*7.5CM	24.5*7.5CM	23*7.5CM
Weight	80G	77G	80G	96G
Material	PVC	PVC	PVC	PVC

Item Number	EYE-GEL 0414A	EYE-GEL 0414B	EYE-GEL 0208	EYE-GB 0208
Description	36*9.5CM	36*9.5CM	21*6CM	21*6CM
Weight	152G	152G	42G	35G
Material	PVA	PVA	PVA	PVC

SPA Series Face Mask

Conforms to your body to deliver soothing relief from:
migraines, sinus infections, cosmetic procedures and more
Freeze or microwave
Reusable and non-toxic
Latex, lead and BPA-free




"Age can always be a Secret"



Item Number	FACE-GEL 0812	FACE-GEL 0812	FACE-GB 0812	FACE-GB 0812
Description	19*30.8CM	19*30.8CM	19*30.8CM	19*30.8CM
Weight	230G	230G	230G	230G
Material	PVC	PVC	PVC	PVC



Disposable Warm Patch Normal Heat



Item Number: HP-BACK
Description: 10*13CM
Weight: 40G
Material: non-woven



Strictly temperature control to ensure safe use
Adjust temperature to meet specific needs
Design for different body parts, enjoy comfortable warm
Disposable and easy to use when doing outdoor activities

Disposable Warm Patch

Using heat patch is the most convenient and safest way to provide skin-touching continuous and stable warm. It is portable and has great effect in preventing from frostbite in winter. Just one single piece can resist full day's cold.

				
Item Number	WP-MUL	WP-HAND	WP-FEET	WP-TOE
Description	10*13CM	9*5.5CM	8*22CM	7*10CM
Weight	40G	26G*2	56G	16G*2
Material	Non-woven	Non-woven	Non-woven	Non-woven



an INTElligent COmpany to partner with

Disposable Warm Patch Pain Relief



Item Number: HP-235
Description: 15*10.5CM
Weight: 23.5G
Material: non-woven



Main features:

It's disposable and easy to provide relief of minor muscular and joint aches and pains associated with over exertion, strains, sprains and arthritis, period pain, etc. It will gradually get warm and last for at least 8 hours.

Instruction:

Open the outer pouch and shake gently, leaving the internal pack to react with air. It needs around 30 minutes to reach its therapeutic temperature, use it on the affected area.

Item Number	HP 25B	HP 82	HP 50	HP 40B/HP 12
Description	9.5*13CM	29.5*9.5CM	29.5*9CM	12*18CM/7*9.5CM
Weight	25G	82G	50G	40G/12G
Material	non-woven	non-woven	non-woven	non-woven
Using Position	For Hypogastrium	For Back and Shoulder	For Back and Shoulder	For Hypogastrium

Item Number	HP 25A	HP 38	WP 186	WP 20
Description	9.5*13CM	22*9.5CM	15*10.5CM	11*9.5CM
Weight	25G	38G	18.6G	20G
Material	non-woven	non-woven	non-woven	non-woven
Using Position	For Upper-Abdomen	For Knee	For Invisible Knee	For Pelma



Home Care Series
More Convenient Life



Home Care Series ICE Bag



Characteristic:

Excellent cold therapy effects.
Pleated design allows bag fabric to conform closely to body contours.
Different sizes and colors available.
Portable and reusable.

Usage:

Use with ice cubes, crushed ice or cold water.
Inject the ice or cold water into the bag and close the bag tightly, then apply it on the body.

Item Number	ICE-BAG 6	ICE-BAG 9	ICE-BAG 11
Description	6 inch	9 inch	11 inch
Weight	68G	97G	120G
Material	polyester	polyester	polyester



Home Care Series Fever Cooling Patch



Item Number: FCP 0305
Description: 5*12CM
Weight: 12.5G
Material: Non-woven



Soothe fever either for Baby or Adult
Relief pain caused by headache or toothache
Keep cool when driving or doing outdoor activities in summer

Advantages:

Gentle to skin , Strong Adhesion , Disposable , Menthol Aroma
Contain Natural non-drug ingredients , in hydro-gel , Safe to use with medications.
Soothe fever instantly , Cooling effect can last up to 8 hours




Reusable Cooling Pack

Reusable , keep in freezer for less than 2 hours before use

Item Number	CP 120	CP 270	CP
Description	10*15CM	13*20CM	13*20CM
Weight	120G	270G	280G
Material	Polyester taffeta/Fabric	EVA	PVC compound spunlace

Home Care Series Cool Gel Mat



Instruction:

For pillow pad, just put it on your pillow. For bed mat, you can lie on it directly or put the bed sheet on the gel mat, then lie on it.

Helps relieve night sweats and hot flashes
 Non-toxic and environmentally friendly
 Adopts special techniques, non-leaking and safe
 Enjoy coolness and softness at the same time
 Keep cool automatically, no electricity required, economical

Item Number	IBD 3040A	IBD 3040B	IBD 4050
Description	30*40CM	30*40CM	40*50CM
Weight	680G	680G	1080G
Material	Cotton/Polyester taffeta/Fabric	Cotton/Polyester taffeta/Fabric	Cotton/Polyester taffeta/Fabric

Item Number	IBD 2030A	IBD 2030B	IBD 3040C
Description	20*30CM	20*30CM	30*40CM
Weight	340G	320G	680G
Material	Cotton/Polyester taffeta/Fabric	Cotton/Polyester taffeta/Fabric	Cotton/Polyester taffeta/Fabric

Item Number	IBD 6090	IBD 9090	IBD 90140
Description	60*90CM	90*90CM	90*140CM
Weight	2760G	4100G	6500G
Material	Cotton/Polyester taffeta/Fabric	Cotton/Polyester taffeta/Fabric	Cotton/Polyester taffeta/Fabric

Home Care Series PCM Cooling Mat

This kind of cooling mat is made of constant temperature phase change materials
 At the room temperature about 26 °C, the cooling effect can last almost 6 hours

Item Number	PCM 3035	PCM 3040	PCM 90140
Description	30*35CM	30*40CM	90*140CM
Weight	620G	680G	6620G
Material	Polyester taffeta/Fabric	Polyester taffeta/Fabric	Polyester taffeta/Fabric

Home Care Series Cooling Gel Cushion



Keep you cool while driving in hot days to ensure your safety
Filled with special sponge and gel, makes you comfortable and cool
No Refrigeration or Electricity Required, Environmentally friendly

At home or At the office
Under the temperature around 26 C , the cushion has a better cooling effect
Extra soft prevents your waist from getting stiffness while working, studying or watching TV




			
Item Number	CUS 2030	CUS 3040	CUS 3035
Description	20*30CM	30*40CM	30*35CM
Weight	860G	1565G	1370G
Material	Polyester taffeta+Sponge	Polyester taffeta+Sponge	Polyester taffeta+Sponge




Home Care Series Cooling Gel Pillow

Advantages:

Lower your body temperature immediately
The cooling effect can lasts up to 3 hours
Non-toxic, Environmental-friendly, Energy-saving
Recent sleep studies have shown that cooling the brain is a simple, effective method to help relieve insomnia.
Soothes sore, aching muscles, numbness, headaches and stress pain. Also great for hot flashes and night sweaters
Absorb heat from the body to help you sleep better in hot days
Ideal choice for home beds, toddler beds, hospital settings lounging around



			
Item Number	CP 650	CP 1000	CP 1040A
Description	18*29.5CM	18.5*31CM	21.6*31.8CM
Weight	650G	1000G	1040G
Material	Polyester taffeta	PE&NYLON+EVA	PVC+PE&NYLON

			
Item Number	CP 400	CP 420A	CP 420B
Description	45.7*32.6CM	54.7*31.4CM	47.5*32.5CM
Weight	400G	420G	420G
Material	PVC	PVC	PVC

Home Care Series Cooling Scarf

The inner macromolecule can absorb large amount of water to lower down the body temperature, bring you a cool feeling when you do outdoor activities.

Instruction:

Put the solid part of this product into water around 10 minutes. If inner material flows, use your hand to make it even. Use towel to absorb extra water. Don't screw it heavily, which may damage the product. Put it on forehead, neck or other suitable position. If it's not cool enough, please use other side or put into water again. Do not freeze it.



		
Item Number	SCA 23	SCA 20
Description	5*105CM	5*90CM
Weight	23G	20G
Material	Cotton	Polyester taffeta

Cooling Belt

Reusable, store in the freezer about 1 hour prior to use Keep cool when doing sports, working outside in hot days For soothing fever, relieve pain cause by headache

		
Item Number	BELT 280	BELT 190
Description	13*46CM	8*55CM
Weight	280G	190G
Material	PE&NYLON	PE&NYLON

Home Care Series PVA Towel



With good hygroscopic property, it is extra soft after absorb water, so it can be perfectly used as:
Cooling Towel to keep cool when doing sports or working outdoors
Cleaning Towel to clean cars, windows, kitchen-ware



		
Item Number	PVA-TW	PVA-TWL
Description	All kinds of sizes	all kinds of sizes
Weight		
Material	PVA	PVA

Mom & Baby Series Breast Gel Pack



Advantages:

- Hot and cold therapy when Breast-feeding or during menstruation
- Cold therapy to relieve pain, swelling, engorgement, plugged ducts, or mastitis
- Hot therapy to help relieve mastitis and help loosen clogged milk ducts to make feeding easier
- Breast Pump Aid: use hot to reduce time spent pumping


Breast Gel Pack

Non-toxic, phthalate free, BPA free, latex free
Soft and flexible for comfort and feature a contoured shape that molds to your breast
With extra soft covers conform to breast and can be worn under bra.


		
Item Number	BRA-GEL	BRA-GB
Description	17*17CM	17*17CM
Weight	165G	165G
Material	PVA	PVA



Heat in water



Heat in microwave oven



Freeze in refrigerator

Mom & Baby Series Perineal Instant Cold Pack

Advantages:

Perineal Instant Cold Pack especially used for gentle, effective postpartum care
Maximum comfort and protection with minimal nursing intervention

Thermal therapy and absorbent pad in one
The all-in-one design eliminates the need for additional supplies, offering superior cost efficiency.

Made of soft material, our perineal packs provide safe, controlled and effective thermal therapy and increased absorbency



Item Number: PER 0613
Description: 14.5*33.5CM
Weight: 104G
Material: PE

Alternate Activation:

1. Position hands on inner water bag.
2. Firmly but gently twist inner bag to rupture.
3. Shake to mix contents.
4. Remove outer packaging and peel back adhesive liner.
5. Apply non-adhesive side to patient and secure pack to undergarment.

Mom & Baby Series Infant Heel Warmers

Heel warmers increase blood flow to infant's heel area, minimizing potential bruising, infection or cartilage damage caused by frequent heel sticks from blood tests.
Specially designed gel pack conforms gently to infant's heel-activates instantly and has lift-and-stick tape for easy application and removal.

Advantages:

Strictly temperature control , gentle to baby skin
Disposable, Safe and no cross-infection



Item Number: IHW 0305
Description: 12*8.5CM
Weight: 57G
Material: PE

Instruction:

Bend metal disk back and forth with thumb and forefinger of both hands. Knead pack and secure warmer to heel by wrapping strip around foot.

Promotion & Gift Series

It's very convenient to use hand warmer outdoors, such as hiking. Bring you warm in the cold weather. It can be made into various shape with different logo to meet your special need. It's good choice as gifts and promotional items with cheap price.

Warm Toy

				
Item Number	HC-GEL 40	HC-GEL 80	HC-GEL 50	HC-GEL 85
Description	9.4*9CM	10.7*10.7CM	9.5*10CM	11.5*11.5CM
Weight	40G	80G	50G	85G
Material	PVC	PVC	PVC	PVC

				
Item Number	HC-GEL350	MAG 70A	MAG 100A	MAG 100B
Description	15.3*22.3CM	6.7*12CM	9.8*9.9CM	10*10CM
Weight	350G	70G	100G	100G
Material	PE&NYLON	PVC	PVC	PVC

Magic Warmer

Instruction:

Bend the metal disk inside the pack, it will start crystallization. When the crystallization is completed, it's ready to use. To recharge, please place hand warmer in boiling water with cloth. Heat it about 5 minutes or until the crystal is dissolved. It's reusable.

		
Item Number	MAG-HEART	MAG-STAR
Description	10.6*11.5CM	9.5*10CM
Weight	100G	53G
Material	PVC	PVC

